

# Starling ACCESSORIES CATALOG



# Table of Contents

1.	Starling	1
2.	Starling Sensors	2
3.	Starling Patient Cable	3
4.	Starling Non-Invasive Blood Pressure Module & Accessories	4
<b>5</b> .	Starling Sp02 Module & Accessories	5
6.	Passive Leg Raise (PLR) LIFT	6
<b>7.</b>	Starling Rapid Response Team Bag	7
8.	Starling Pole Mount	8
9.	Starling Roll Stand	9
10.	Starling Simulator	10
11.	Connectivity Accessories	10
12.	Starling Battery Pack	11
13.	Starling Power Cable	11

# Starling

### Non-Invasive Fluid Management Monitoring System

Starling Monitor
\*includes Starling Patient Cable & Starling Power Cable

P/N: CMMST5



The **Starling** non-invasive fluid management monitoring system provides continuous hemodynamic monitoring and empowers fluid management across the continuum of care: **ICU, ED, RRT and OR.** 

A portable and lightweight touch screen monitor:

**H:** 8.7", **W:** 11.4", **D:** 7.4", 9.5 lb, 10.4"

6 Hr. Battery Life

# Starling Sensors

Sensors — Box of 10 P/N: CMS10

Sensors — Box of 25 P/N: CMS25

Sensors — Box of 50 P/N: CMS50



Starling sensors can be placed anywhere on the chest or back as long as two are positioned above and below the heart.
Sensors should be replaced every 48 hours.



# 3 Starling Patient Cable

Patient Cable

P/N: CMASC10
\*one included with **Starling** Monitor



The **Starling** patient cable connects the **Starling** sensors to the **Starling** monitor.

Length:  $310 \pm 7$  cm

Material: Covered with PVC

### Starling Non-Invasive Blood Pressure (NIBP) Module & Accessories

Starling Blood Pressure Module
\*includes Standard Adult Cuff and Hose

P/N: CMANIBP

#### Accessories

NIBP Hose

P/N: CMAHOSE
\*one included with module

Standard Adult NIBP Cuff

P/N: CMACU23
\*one included with module

Large Adult NIBP Cuff

P/N: CMACU31



The **Starling** blood pressure module enables non-invasive blood pressure monitoring on the **Starling** monitor. The **Starling** blood pressure module includes a reusable SunTech Medical standard adult NIBP cuff and hose. Additional cuff and hose accessories can also be ordered separately.

# 5 Starling Sp02 Module & Accessories

**Starling** Sp02 Module \*includes Cable and Reusable Finger Probe Set

P/N: CMASP02

#### **Accessories**

Cable and Reusable Finger Probe Set

P/N: CMANELLEP \*one included with module



The **Starling** Sp02 module enables Sp02 monitoring on the **Starling** monitor. The **Starling** Sp02 module includes a Medtronic Sp02 cable and reusable finger probe set. Additional finger probe sets can be ordered separately.

# 6 Passive Leg Raise (PLR) LIFT

LIFT

P/N: CMALIFT



Portable passive leg raise lift, opens at a ~45 degree angle

**Size:** 50x60x80 cm

# Starling Rapid Response Team Bag

Starling RRT Bag P/N: CMARRTBAG




The **Starling** RRT bag is a portable bag that holds all the essentials for treating a hypotension emergency, including: Sp02, BP and IV access supplies.

The **Starling** monitor can be operated through the transparent front cover.

# Starling Pole Mount

Starling Pole Mount

P/N: CMASPM30



Enables mounting the **Starling** monitor to an IV pole.

### Starling Roll Stand

Starling Roll Stand

P/N: CMASRS20S1



**Starling** roll stand (GCX) provides optimal, secure and convenient positioning of the **Starling** monitor within the clinical environment. Equipped with a handle for easy maneuverability and five caster wheels with lockable brakes.

# Starling Simulator

Simulator

P/N: CMASSIM



### 11 Connectivity Accessories

Philips connectivity Cables and connector set
Used to connect the **Starling** monitor
with Philips Patient Monitor

GE connectivity Cables and connector set
Used to connect the **Starling** monitor with
GE Patient Monitor.

# 12 Starling Battery Pack 13 Starling Power Cable

Replacement Battery Pack P/N: CMABATT

Power Cable

P/N: CPA1705 \*one included with **Starling** Monitor

5m length power cable





**Rx Only.** For safe and proper use of product mentioned herein, please refer to the Instructions for Use or Operators Manual.

#### Baxter.com

Baxter International inc. One Baxter Parkway / Deerfield, Illinois 60015

Baxter and Starling are trademarks of Baxter International Inc. or its subsidiaries.

USMP/CHE/20-0013 05/20