

Baxter

Starling

FLUID MANAGEMENT
MONITORING SYSTEM

Starling

ACCESSORIES CATALOG

Table of Contents

1. Starling	1
2. Starling Sensors	2
3. Starling Patient Cable	3
4. Starling Non-Invasive Blood Pressure Module & Accessories	4
5. Starling SpO2 Module & Accessories	5
6. Passive Leg Raise (PLR) LIFT	6
7. Starling Rapid Response Team Bag	7
8. Starling Pole Mount	8
9. Starling Roll Stand	9
10. Starling Simulator	10
11. Connectivity Accessories	10
12. Starling Battery Pack	11
13. Starling Power Cable	11

1 Starling

Non-Invasive Fluid Management Monitoring System

Starling Monitor

*Includes **Starling** Patient Cable & **Starling** Power Cable

P/N: CMMST5



The **Starling** non-invasive fluid management monitoring system provides continuous hemodynamic monitoring and empowers fluid management across the continuum of care: **ICU, ED, RRT and OR.**

A portable and lightweight touch screen monitor:

H: 8.7", W: 11.4", D: 7.4",
9.5 lb, 10.4"

6 Hr. Battery Life

2

Starling Sensors

Sensors — Box of 10 P/N: CMS10

Sensors — Box of 25 P/N: CMS25



Starling sensors can be placed anywhere on the chest or back as long as two are positioned above and below the heart. Sensors should be replaced every 48 hours.



3

Starling Patient Cable

Patient Cable

P/N: CMASC10

*one included with **Starling** Monitor



The **Starling** patient cable connects the **Starling** sensors to the **Starling** monitor.

Length: 310 ± 7 cm

Material: Covered with PVC

4 **Starling** Non-Invasive Blood Pressure [NIBP] Module & Accessories

Starling Blood Pressure Module

*includes Standard Adult Cuff and Hose

P/N: CMANIBP

Accessories

NIBP Hose

P/N: CMAHOSE

*one included with module

Standard Adult NIBP Cuff

P/N: CMACU23

*one included with module

Large Adult NIBP Cuff

P/N: CMACU31



The **Starling** blood pressure module enables non-invasive blood pressure monitoring on the **Starling** monitor. The **Starling** blood pressure module includes a reusable SunTech Medical standard adult NIBP cuff and hose. Additional cuff and hose accessories can also be ordered separately.

5 **Starling** SpO2 Module & Accessories

Starling SpO2 Module

*includes Cable and Reusable Finger Probe Set

P/N: CMASPO2

Accessories

Cable and Reusable Finger
Probe Set

P/N: CMANELLEP

*one included with module



The **Starling** SpO2 module enables SpO2 monitoring on the **Starling** monitor. The **Starling** SpO2 module includes a Medtronic SpO2 cable and reusable finger probe set. Additional finger probe sets can be ordered separately.

6 Passive Leg Raise (PLR) LIFT

LIFT P/N: CMALIFT



Portable passive leg raise lift,
opens at a ~45 degree angle

Size: 50x60x80 cm

7

Starling Rapid Response Team Bag

Starling RRT Bag P/N: CMARRTBAG




The **Starling** RRT bag is a portable bag that holds all the essentials for treating a hypotension emergency, including: SpO2, BP and IV access supplies.

The **Starling** monitor can be operated through the transparent front cover.

8

Starling Pole Mount

Starling Pole Mount P/N: CMASPM30



Enables mounting the **Starling** monitor to an IV pole.

9

Starling Roll Stand

Starling Roll Stand P/N: CMASRS20S1



Starling roll stand (GCX) provides optimal, secure and convenient positioning of the **Starling** monitor within the clinical environment. Equipped with a handle for easy maneuverability and five caster wheels with lockable brakes.

10 **Starling** Simulator

Simulator P/N: CMASSIM



11 Connectivity Accessories

Philips connectivity Cables and connector set P/N: CCA3601PH


Used to connect the **Starling** monitor with Philips Patient Monitor

GE connectivity Cables and connector set P/N: CCA3601GE

Used to connect the **Starling** monitor with GE Patient Monitor.

12 **Starling** Battery Pack

Replacement Battery Pack P/N: CMABATT



13 **Starling** Power Cable

Power Cable P/N: CPA1705
*one included with **Starling** Monitor

5m length power cable



Rx Only. For safe and proper use of product mentioned herein,
please refer to the Instructions for Use or Operators Manual.

Baxter.com

Baxter International inc.
One Baxter Parkway / Deerfield, Illinois 60015

Baxter and Starling are trademarks of
Baxter International Inc. or its subsidiaries.

USMP/CHE/20-0013 05/20